

TBR2 Polyclonal Antibody Bio Conjugated

Catalog No: #C90677Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TBR2 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	TBR2
Other Names	TBR2
Accession No.	Uniprot:O95936Gene ID:8320
Uniprot	O95936
GeneID	8320
Calculated MW	73
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only