

CLSTN2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90680Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLSTN2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	261-360/955
Conjugates	Bio
Target Name	CLSTN2
Other Names	CLSTN2
Accession No.	Uniprot:Q9H4D0Gene ID:64084
Uniprot	Q9H4D0
GeneID	64084
Calculated MW	105
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only