

LAR Polyclonal Antibody PE Conjugated

Catalog No: #C90708PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LAR Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	30-100/1907
Conjugates	PE
Target Name	LAR
Other Names	LAR
Accession No.	Uniprot:P10586Gene ID:5792
Uniprot	P10586
GeneID	5792
Calculated MW	210
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only