

## Fibrillin 1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90715Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fibrillin 1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/2872
Conjugates	Bio
Target Name	Fibrillin 1
Other Names	Fibrillin 1
Accession No.	Uniprot:P35555Gene ID:2200
Uniprot	P35555
GeneID	2200
Calculated MW	316
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only