

AVPR1A Polyclonal Antibody Cy3 Conjugated

Catalog No: #C90718Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AVPR1A Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-55/418
Conjugates	Cy3
Target Name	AVPR1A
Other Names	AVPR1A
Accession No.	Uniprot:P37288Gene ID:552
Uniprot	P37288
GeneID	552
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only