

# Nuclear Pore Complex Proteins/Nup107 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90733Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Nuclear Pore Complex Proteins/Nup107 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	501-600/925
Conjugates	Bio
Target Name	Nuclear Pore Complex Proteins/Nup107
Other Names	Nuclear Pore Complex Proteins/Nup107
Accession No.	Uniprot:P57740Gene ID:57122
Uniprot	P57740
GeneID	57122
Calculated MW	106
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only