

NPEPPS Polyclonal Antibody PE Conjugated

Catalog No: #C90735PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------------------------|
| Product Name | NPEPPS Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 201-300/919 |
| Conjugates | PE |
| Target Name | NPEPPS |
| Other Names | NPEPPS |
| Accession No. | Uniprot:P55786 Gene ID:9520 |
| Uniprot | P55786 |
| GeneID | 9520 |
| Calculated MW | 103 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only