

Tetanus Toxin heavy chain Polyclonal Antibody FITC Conjugated

Catalog No: #C90749FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tetanus Toxin heavy chain Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Clostridium tetani E88
Immunogen Type	801-900/1315
Conjugates	FITC
Target Name	Tetanus Toxin heavy chain
Other Names	Tetanus Toxin heavy chain
Accession No.	Uniprot:P04958Gene ID:24255210
Uniprot	P04958
GeneID	24255210
Calculated MW	94
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only