

PAR-2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90753Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAR-2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-150/397
Conjugates	Bio
Target Name	PAR-2
Other Names	PAR-2
Accession No.	Uniprot:P55085Gene ID:2150
Uniprot	P55085
GeneID	2150
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only