

PTF1A Polyclonal Antibody Biotin Conjugated

Catalog No: #C90775Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | PTF1A Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IHC ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 165-260/328 |
| Conjugates | Bio |
| Target Name | PTF1A |
| Other Names | PTF1A |
| Accession No. | Uniprot:Q7RTS3Gene ID:256297 |
| Uniprot | Q7RTS3 |
| GeneID | 256297 |
| Calculated MW | 35 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only