

# HLA-DR/HLA DRB1 Polyclonal Antibody PE-Cy5 Conjugated



Catalog No: #C90783PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HLA-DR/HLA DRB1 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/254
Conjugates	PE-Cy5
Target Name	HLA-DR/HLA DRB1
Other Names	HLA-DR/HLA DRB1
Accession No.	Uniprot:P01903Gene ID:3122
Uniprot	P01903
GeneID	3122
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only