

PAX8 Polyclonal Antibody FITC Conjugated

Catalog No: #C90789FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PAX8 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC IF
Species Reactivity	Hu Ms Rat
Conjugates	FITC
Target Name	PAX8
Other Names	PAX8
Accession No.	Uniprot:Q06710Gene ID:7849
Uniprot	Q06710
GeneID	7849
Calculated MW	48
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only