

GPRC5B Polyclonal Antibody PE Conjugated

Catalog No: #C90791PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPRC5B Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/403
Conjugates	PE
Target Name	GPRC5B
Other Names	GPRC5B
Accession No.	Uniprot:Q9NZH0Gene ID:51704
Uniprot	Q9NZH0
GeneID	51704
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only