

## GLI1 Polyclonal Antibody PE Conjugated

Catalog No: #C90802PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GLI1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/1106
Conjugates	PE
Target Name	GLI1
Other Names	GLI1
Accession No.	Uniprot:P08151Gene ID:2735
Uniprot	P08151
GeneID	2735
Calculated MW	118
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only