

P2Y12 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90805Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2Y12 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	141-240/342
Conjugates	Bio
Target Name	P2Y12
Other Names	P2Y12
Accession No.	Uniprot:Q9H244Gene ID:64805
Uniprot	Q9H244
GeneID	64805
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only