

ACAN Polyclonal Antibody APC Conjugated

Catalog No: #C90827APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ACAN Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-220/2415
Conjugates	APC
Target Name	ACAN
Other Names	ACAN
Accession No.	Uniprot:P16112Gene ID:176
Uniprot	P16112
GeneID	176
Calculated MW	208/248
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only