

Adenosine A3 Receptor Polyclonal Antibody PE Conjugated

Catalog No: #C90829PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adenosine A3 Receptor Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC IF
Species Reactivity	Hu Ms Rat
Conjugates	PE
Target Name	Adenosine A3 Receptor
Other Names	Adenosine A3 Receptor
Accession No.	Uniprot:P0DMS8Gene ID:140
Uniprot	P0DMS8
GeneID	140
Calculated MW	35
Concentration	1 mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only