

# Phospho-TIE2 (Tyr1108) Polyclonal Antibody Biotin Conjugated



Catalog No: #C90830Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospho-TIE2 (Tyr1108) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TL(p-Y)EK
Conjugates	Bio
Target Name	Phospho-TIE2 (Tyr1108)
Other Names	Phospho-TIE2 (Tyr1108)
Accession No.	Uniprot:Q02763Gene ID:7010
Uniprot	Q02763
GeneID	7010
Calculated MW	124
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only