

Fragilis Polyclonal Antibody PE Conjugated

Catalog No: #C90831PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fragilis Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	FC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	41-100/133
Conjugates	PE
Target Name	Fragilis
Other Names	Fragilis
Accession No.	Uniprot:Q01628 Gene ID:10410
Uniprot	Q01628
GeneID	10410
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only