

ABCG1 Polyclonal Antibody PE Conjugated

Catalog No: #C90837PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ABCG1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	561-678/678
Conjugates	PE
Target Name	ABCG1
Other Names	ABCG1
Accession No.	Uniprot:P45844Gene ID:9619
Uniprot	P45844
GeneID	9619
Calculated MW	76
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only