

Stabilin 2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90841Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Stabilin 2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2201-2300/2551
Conjugates	Bio
Target Name	Stabilin 2
Other Names	Stabilin 2
Accession No.	Uniprot:Q8WWQ8Gene ID:55576
Uniprot	Q8WWQ8
GeneID	55576
Calculated MW	154/275
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only