

PNPLA3 Polyclonal Antibody PE Conjugated

Catalog No: #C90878PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PNPLA3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu
Immunogen Type	101-200/481
Conjugates	PE
Target Name	PNPLA3
Other Names	PNPLA3
Accession No.	Uniprot:Q9NST1Gene ID:80339
Uniprot	Q9NST1
GeneID	80339
Calculated MW	53
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only