

ZNF385B Polyclonal Antibody FITC Conjugated

Catalog No: #C90885FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZNF385B Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-100/471
Conjugates	FITC
Target Name	ZNF385B
Other Names	ZNF385B
Accession No.	Uniprot:Q569K4Gene ID:151126
Uniprot	Q569K4
GeneID	151126
Calculated MW	50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only