

phospho-beta 2 Adrenergic Receptor (Ser355 / Ser356) Polyclonal Antibody Biotin Conjugated

Catalog No: #C90902Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-beta 2 Adrenergic Receptor (Ser355 / Ser356) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	GY(p-S)(p-S)N
Conjugates	Bio
Target Name	phospho-beta 2 Adrenergic Receptor (Ser355 / Ser356)
Other Names	phospho-beta 2 Adrenergic Receptor (Ser355 / Ser356)
Accession No.	Uniprot:P07550Gene ID:154
Uniprot	P07550
GeneID	154
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only