

## Syndecan-1 Polyclonal Antibody APC Conjugated

Catalog No: #C90912APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Syndecan-1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-310/310
Conjugates	APC
Target Name	Syndecan-1
Other Names	Syndecan-1
Accession No.	Uniprot:P18827Gene ID:6382
Uniprot	P18827
GeneID	6382
Calculated MW	34
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only