

# Integrin Alpha V / Beta 3 (CD51/CD61) Polyclonal Antibody Biotin Conjugated

Catalog No: #C90914Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin Alpha V / Beta 3 (CD51/CD61) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	Integrin Alpha V / Beta 3 (CD51/CD61)
Other Names	Integrin Alpha V / Beta 3 (CD51/CD61)
Accession No.	Uniprot:P05106, P06756Gene ID:3685
Uniprot	P05106, P06756
GeneID	3685
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only