

LYVE-1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90915Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	LYVE-1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	231-318/318
Conjugates	Bio
Target Name	LYVE-1
Other Names	LYVE-1
Accession No.	Uniprot:Q9Y5Y7Gene ID:114332
Uniprot	Q9Y5Y7
GeneID	114332
Calculated MW	32
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only