

phospho-FANCD2 (Ser1404) Polyclonal Antibody PE Conjugated

Catalog No: #C90922PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-FANCD2 (Ser1404) Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu
Immunogen Type	QN(p-S)QE
Conjugates	PE
Target Name	phospho-FANCD2 (Ser1404)
Other Names	phospho-FANCD2 (Ser1404)
Accession No.	Uniprot:Q9BXW9Gene ID:2177
Uniprot	Q9BXW9
GeneID	2177
Calculated MW	162
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only