

GAA Polyclonal Antibody Biotin Conjugated

Catalog No: #C90938Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | GAA Polyclonal Antibody Biotin Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IHC ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 541-640/952 |
| Conjugates | Bio |
| Target Name | GAA |
| Other Names | GAA |
| Accession No. | Uniprot:P10253Gene ID:2548 |
| Uniprot | P10253 |
| GeneID | 2548 |
| Calculated MW | 83/92/98 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only