

Metallothionein Polyclonal Antibody AP Conjugated

Catalog No: #C90943AP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Metallothionein Polyclonal Antibody AP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IHC |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 11-61/61 |
| Conjugates | AP |
| Target Name | Metallothionein |
| Other Names | Metallothionein |
| Accession No. | Uniprot:Q63871Gene ID:17749 |
| Uniprot | Q63871 |
| GeneID | 17749 |
| Calculated MW | 6 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only