

ZO-1 Polyclonal Antibody APC Conjugated

Catalog No: #C90946APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ZO-1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1551-1702/1733
Conjugates	APC
Target Name	ZO-1
Other Names	ZO-1
Accession No.	Uniprot:Q07157Gene ID:7082
Uniprot	Q07157
GeneID	7082
Calculated MW	191
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only