

# phospho-GFAP (Ser38) Polyclonal Antibody Biotin Conjugated



Catalog No: #C90948Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-GFAP (Ser38) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms
Immunogen Type	RL(p-S)L
Conjugates	Bio
Target Name	phospho-GFAP (Ser38)
Other Names	phospho-GFAP (Ser38)
Accession No.	Uniprot:P14136Gene ID:2670
Uniprot	P14136
GeneID	2670
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only