

Beta galactosidase Polyclonal Antibody FITC Conjugated

Catalog No: #C90951FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta galactosidase Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/677
Conjugates	FITC
Target Name	Beta galactosidase
Other Names	Beta galactosidase
Accession No.	Uniprot:P16278Gene ID:2720
Uniprot	P16278
GeneID	2720
Calculated MW	73
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=3ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only