

GLEPP1 Polyclonal Antibody Bio Conjugated

Catalog No: #C90952Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLEPP1 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	GLEPP1
Other Names	GLEPP1
Accession No.	Uniprot:Q16827Gene ID:5800
Uniprot	Q16827
GeneID	5800
Calculated MW	135
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC 1:100-500

IF 1:100-500

Note: This product is for in vitro research use only