

FAK Polyclonal Antibody FITC Conjugated

Catalog No: #C90959FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FAK Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	901-1052/1052
Conjugates	FITC
Target Name	FAK
Other Names	FAK
Accession No.	Uniprot:Q05397Gene ID:5747
Uniprot	Q05397
GeneID	5747
Calculated MW	116
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only