

Fibronectin/Ugl-Y3 Polyclonal Antibody FITC Conjugated

Catalog No: #C90963FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fibronectin/Ugl-Y3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	2201-2300/2386
Conjugates	FITC
Target Name	Fibronectin/Ugl-Y3
Other Names	Fibronectin/Ugl-Y3
Accession No.	Uniprot:P02751Gene ID:2335
Uniprot	P02751
GeneID	2335
Calculated MW	259
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only