

## Granzyme B Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90970PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Granzyme B Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-130/247
Conjugates	PE-Cy5
Target Name	Granzyme B
Other Names	Granzyme B
Accession No.	Uniprot:P10144Gene ID:3002
Uniprot	P10144
GeneID	3002
Calculated MW	25
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only