

MCL1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90971Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MCL1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-350/350
Conjugates	Bio
Target Name	MCL1
Other Names	MCL1
Accession No.	Uniprot:Q07820Gene ID:4170
Uniprot	Q07820
GeneID	4170
Calculated MW	39
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only