

GPR43 Polyclonal Antibody Bio Conjugated

Catalog No: #C90975Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPR43 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	GPR43
Other Names	GPR43
Accession No.	Uniprot:O15552Gene ID:2867
Uniprot	O15552
GeneID	2867
Calculated MW	37KD
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only