

CLEC9A Polyclonal Antibody APC Conjugated

Catalog No: #C90985APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CLEC9A Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Rat
Immunogen Type	101-200/241
Conjugates	APC
Target Name	CLEC9A
Other Names	CLEC9A
Accession No.	Uniprot:Q6UXN8Gene ID:283420
Uniprot	Q6UXN8
GeneID	283420
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only