

## TMEM176A Polyclonal Antibody PE Conjugated

Catalog No: #C90986PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TMEM176A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms
Immunogen Type	1-100/235
Conjugates	PE
Target Name	TMEM176A
Other Names	TMEM176A
Accession No.	Uniprot:Q96HP8Gene ID:55365
Uniprot	Q96HP8
GeneID	55365
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only