

AIF1/Iba1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90987Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AIF1/Iba1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-147/147
Conjugates	Bio
Target Name	AIF1/Iba1
Other Names	AIF1/Iba1
Accession No.	Uniprot:P55008Gene ID:199
Uniprot	P55008
GeneID	199
Calculated MW	16
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only