

SLAP Polyclonal Antibody PE Conjugated

Catalog No: #C90988PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLAP Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/276
Conjugates	PE
Target Name	SLAP
Other Names	SLAP
Accession No.	Uniprot:Q13239Gene ID:6503
Uniprot	Q13239
GeneID	6503
Calculated MW	31
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only