

SSTR3 Polyclonal Antibody FITC Conjugated

Catalog No: #C90989FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SSTR3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	121-230/418
Conjugates	FITC
Target Name	SSTR3
Other Names	SSTR3
Accession No.	Uniprot:P32745Gene ID:6753
Uniprot	P32745
GeneID	6753
Calculated MW	46
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only