

# Centrin-1/Centrin-2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C91007Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Centrin-1/Centrin-2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-172/172
Conjugates	Bio
Target Name	Centrin-1/Centrin-2
Other Names	Centrin-1/Centrin-2
Accession No.	Uniprot:Q12798Gene ID:1068 1069
Uniprot	Q12798
GeneID	1068 1069
Calculated MW	20
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only