

ERPINA6/Cortisol Binding Globulin Polyclonal Antibody HRP Conjugated



Catalog No: #C91019HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ERPINA6/Cortisol Binding Globulin Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms
Immunogen Type	121-220/405
Conjugates	HRP
Target Name	ERPINA6/Cortisol Binding Globulin
Other Names	ERPINA6/Cortisol Binding Globulin
Accession No.	Uniprot:P08185Gene ID:866
Uniprot	P08185
GeneID	866
Calculated MW	43
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only