

CSAG1 Polyclonal Antibody HRP Conjugated

Catalog No: #C91026HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CSAG1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu
Immunogen Type	22-78/78
Conjugates	HRP
Target Name	CSAG1
Other Names	CSAG1
Accession No.	Uniprot:Q6PB30Gene ID:158511
Uniprot	Q6PB30
GeneID	158511
Calculated MW	6.7
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only