

Histone H1t Polyclonal Antibody Biotin Conjugated

Catalog No: #C91032Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone H1t Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	51-207/207
Conjugates	Bio
Target Name	Histone H1t
Other Names	Histone H1t
Accession No.	Uniprot:P22492Gene ID:3010
Uniprot	P22492
GeneID	3010
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200

Note: This product is for in vitro research use only