

AHR Polyclonal Antibody APC Conjugated

Catalog No: #C91034APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AHR Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	10-120/948
Conjugates	APC
Target Name	AHR
Other Names	AHR
Accession No.	Uniprot:P35869Gene ID:196
Uniprot	P35869
GeneID	196
Calculated MW	93
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only