

EPOR Polyclonal Antibody PE Conjugated

Catalog No: #C91037PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EPOR Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-450/508
Conjugates	PE
Target Name	EPOR
Other Names	EPOR
Accession No.	Uniprot:P19235Gene ID:2057
Uniprot	P19235
GeneID	2057
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only